



A MOBILE WEIGHING CART FOR ROOMING-IN

CUSTOMIZABLE | SAFE | MOBILE

Scale carts are used when rooming-in is strongly practiced; in some cases, the nursery has been repurposed, requiring in-room weighing of newborns. Having a scale cart that looks more like the room's furniture is more appealing to the newborn's family than the cold, clinical stainless steel of traditional scale carts.



BENEFITS

- Designed and built for the 24/7/365 demands of the healthcare environment
- Premium heavy-duty drawer slides with detents to prevent opening while cart is in motion
- Tested to withstand healthcare germicidal solutions
- Standard 5-Year Warranty
- Bumper strip around cart base
- Choice of finishes

STANDARD FEATURES

- Premium precision 5" diameter casters for easy-rolling
 - 2 Locking, 2 Non-Locking

STANDARD DIMENSIONS & CONFIGURATIONS

- 27"W x 22.25"D x 37.5"H
- 2-Drawers
- Door Only
- Drawer Over Doors
- Doors Over Drawers
- Custom Configurations


AVAILABLE OPTIONS

- Thermofoil finish tops and drawer fronts to enhance infection control
- Workflow & safety enhancing options:
 - Single or Dual Flip-up side surfaces
 - Locking drawers/doors: Keyed, push button or electronic lock
 - Drawer dividers
- Custom: Design, Finishes and/or Hardware to accommodate your specific...
 - Workflow: to minimize/eliminate 'workarounds' to enhance safety, accuracy and efficiency
 - Room design: dimensions, finishes, hardware, etc.

Cartstrong Finishes



Kensington Maple



Candlelight



Amber Cherry



Acajou Mahogany



Chocolate Pear Tree